

VaxArray® Influenza Seasonal Hemagglutinin Potency Assay: Array Layout and Specificity


Background: The VaxArray Seasonal Hemagglutinin Potency Assay is a new tool for hemagglutinin protein quantification based on a panel of subtype-specific but broadly reactive monoclonal antibodies (mAbs). Multiple antibodies against seasonal A/H1, A/H3, B/Yamagata-like and B/Victoria-like strains are printed in an array format on a glass substrate. As a multiplexed immunoassay, signal readout is based on fluorescence from a conjugated “universal” primary antibody label.

Objective: The objective of this Technical Note is a basic description of the specificity response to current seasonal strains of influenza virus.

The microarray layout for VaxArray Influenza is illustrated in **Figure 1**. The array contains 9 replicate spots of each monoclonal antibody. There are 3 distinct antibodies for A/H1 and A/H3, 2 antibodies for B/Yamagata-like, and a single antibody for B/Victoria-like strains. The presence of multiple antibodies against the same strain provides breadth of coverage across variations in the virus, as well as complementary insight into protein structure.


Analysis is conducted in a simple sandwich assay format, with a nearly universal, broadly reactive polyclonal antibody (fluor-conjugated) used for detection.

Figure 1 – Array Layout



The specificity of antibodies incorporated in the VaxArray Influenza array toward HA subtypes was evaluated using a range of recombinant HA antigens. **Figure 2** shows representative responses for A/H1, A/H3, B/Yamagata-like, and B/Victoria-like subtypes. Over the range of antigens tested, including recombinant HA from 1999-2011 origin strains, the microarray exhibited excellent specificity and no cross-reactivity between subtypes. It is evident in **Figure 2** that individual H1 strains may respond differently to the panel but, in all tested cases, quantification of HA using at least one antibody was possible.

Figure 2 – Representative Fluorescence Images for Seasonal HA Subtypes and Lineages



The advantage of including multiple distinct mAbs binding to either variable HA1 or conserved HA2 domains can be seen in Figure 2. The capture antibody in position H3 (ii) reliably responded to all of the H3 subtypes originating from viruses first isolated over the time period of 2001 to 2011, despite antigenic drift. Antibody H3 (ii) is known to bind to the conserved HA2 stem region of the protein. The same observation is made for recombinant H1's from viruses isolated over the

time period of 1999 through 2009. The mAb that responded to all H1's (array position H1 (ii)) is known to recognize a conserved region of HA2.

A similar multiple capture approach for each subtype would be impractical for an ELISA system due to reagent cost and low sample throughput per plate (4-6 samples per plate). In contrast, microarrays are ideally suited for the task due to sparing use of reagents (picograms per spot) and multiplex design capability. For example, the amount of antibody needed to coat *a single* ELISA plate can be used to produce over 400 microarrays.

Summary: The VaxArray Influenza assay displays strain specificity, which enables analysis of monovalent or multivalent formulations, as well as robust response despite antigenic drift.