

# VaxArray Assessment of Influenza Vaccine Potency and Stability

Rose Byrne-Nash, Ph.D. | Laura Kuck, Ph.D. | Kathy Rowlen, Ph.D.  
 InDevR, Inc. | 2100 Central Avenue, Suite 106 | Boulder, CO 80301 | Ph : 303.402.9100 | www.indevr.com | E-mail: info@indevr.com

## Abstract




**Background:** There is an on-going effort within the influenza vaccine industry to identify and validate an alternative to the single-radial immunodiffusion assay (SRID) for potency and stability assessment. The VaxArray Seasonal Influenza (VSI) potency assay was developed under the Influenza Vaccine Manufacturing Initiative and performance has been evaluated against key regulatory requirements such as accuracy, precision, and ability to track vaccine stability. This report focuses on a summary of results for a diverse set of flu vaccines tested in collaborative studies.

**Methods:** VSI is based on a panel of monoclonal antibodies printed in a microarray format. A simple sandwich assay is used for simultaneous quantification of all hemagglutinin components within mono- and multi-valent influenza vaccines, including differentiation between the two Influenza B lineages. A standard protocol was used to evaluate the potency of mono-valent and multi-valent influenza vaccines produced in eggs and cell-culture, including recombinant and virus-like particle vaccines. Accuracy was defined with respect to SRID and (or) purity-adjusted total protein content. Potency was also evaluated as a function of time under accelerated stress conditions.






**Conclusions:** VSI was found to exhibit high precision (CV < 15%) and good correlation (R > 0.85) with SRID for most vaccines, and if calibrated appropriately, excellent agreement with SRID values (100 ± 5%). In thermal stability testing all strains in multi-valent vaccines exhibited a decrease in potency over time. While the absolute measured potency for degraded vaccines was not always in agreement with the SRID value, the general trends were consistent with those observed for SRID. These studies indicate that VSI is a reliable and robust alternative potency assay.

## Background

- Potency assays measure the concentration of functional hemagglutinin (HA), which is an influenza virus surface protein. HA has been established to be the key component of whole virus vaccines and the dominant target of protective antibodies following vaccination or infection.
- Currently, the gold standard for influenza vaccine potency is the single radial immunodiffusion (SRID) assay, which has inherent disadvantages including labor-intensive protocols and the requirement for reference reagents that do not necessarily accurately represent the composition of vaccines.
- VaxArray Influenza (VSI) is a rapid alternative potency assay. A comparison between SRID and VSI is demonstrated in Figure 1.



## Results



## The Future of the VaxArray Platform


### Rapid Quantification of Flu Neuraminidase

The VaxArray platform is highly versatile and can be outfitted to screen for a number of different influenza proteins.

Neuraminidase (NA), the second most abundant influenza surface glycoprotein, is gaining support for inclusion in seasonal vaccinations (Eichelberger et al, 2015).

The NA chip is currently under development and will be printed with broadly reactive, subtype specific anti-neuraminidase antibodies (similar to VSI) to A/N1, A/N2, and B-Influenza neuraminidase.

As with VSI, a "universal" polyclonal antibody is used to quantify all neuraminidase in mono-valent and multi-valent mixtures.



Array Position	Detection Target
A	Broad N1
B	Broad N1
C	Broad N1
D	Broad N2
E	Broad N2
F	Broad N2
G	Broad N2
H	Broad B NA
I	Broad B NA

Unlike HA, neuraminidase (NA) does not have a potency gold standard (like SRID for HA), although MUNANA is widely used. We anticipate that this new microarray assay will be able to quantify not only total NA protein content, but also NA enzymatic activity, both important outputs for vaccine production and regulation.

### Application to Pandemic Influenza Vaccines

In concert with the WHO and US government, a number of vaccine producers prepare monovalent vaccines against potentially pandemic flu viruses.

We are developing the VaxArray Pandemic Influenza (VPI) potency assay for A/H5, A/H7, and A/H9 vaccines.

Current studies are ongoing to select for antibodies (shown below) that will ensure that precise, rapid quantification is attainable for high risk pandemic influenza strains.

Antibody ID	Sub-type	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27
1	H5	75% Coverage, 100% Specificity																										

Antibody ID	Sub-type	28	29	30	31	32	33	34	35	36	37	38	
2	H7	64% Coverage, 100% Specificity											
3	H7	64% Coverage, 100% Specificity											

Antibody ID	Subtype	39	40	41	42	43	44	45	46	47	48	49	
4	H9	70% Coverage, 100% Specificity											

VPI will allow rapid, efficient screening of pandemic influenzas as necessary. VPI can function on the leading edge of pandemic vaccine production without requiring reference antisera (similar to VSI), making it a valuable tool for vaccine production, especially during a pandemic. With constant monitoring of potential pandemic influenzas, VPI assays will be adapted to be ready for the next pandemic influenza outbreak.

### Other Advances

- 21 CFR Part 11 capable software is in development to enable easy integration into regulated environments.
- The Vidia Microarray Imaging System is CE certified.
- Compatible with automated and high throughput processing options.
- Please contact us if you are interested in collaborating or in beta testing new products

**Acknowledgements**

We gratefully acknowledge funding from the National Institute for Allergy and Infectious Disease (NIH/NIAID R44 AI102318).

### Characteristics for an Improved Potency Assay for Inactivated Influenza Vaccines

WHO Expert Committee on Biological Standardization; Geneva Oct 2011

Biologically relevant potency measure ✓	Applicability to existing and novel vaccines ✓
Correlation with clinical efficacy (?)	Unaffected by adjuvants (?)
Stability indicating ✓	Measure low doses ✓
Bridging to single radial immunodiffusion (SRD) test ✓	Independent of strain specific reagents X
Precision and accuracy ✓	Reduced amount of reagents (required for assay) ✓
Reproducibility ✓	Robustness of reagent supply, speed of supply, volume/quantity of supply ✓
(Sub)type specificity ✓	Usable in process control (i.e. in presence of other proteins, contaminants) ✓
Flexibility and maximum practicability ✓	Efficient regulatory review ✓
Applicable world-wide ✓	Accelerating lot release ✓
Quick availability/usability following a strain change ✓	