

Good things come
in small
packages



Introducing



ampliPHOX[®]


Colorimetric Detection
Technology



A little gem of technology for microarrays



ampliPHOX[®] provides
fast, visual results




With equivalent sensitivity
to traditional fluorescence



[click here](#)

ampliPHOX[®] uses a light sensitive label that initiates the formation of a polymer where microarray binding events have occurred



1 / 10



of the price of a
fluorescence reader

Results in **15**
minutes



Immediately
visible to the eye



click here

The software images and interprets the final result



The possibilities
are endless

[click here](#)



You develop the
array



ampliPHOX[®]
delivers
the results



Brought to you by



Inspired Life Science Technology

An Environmentally Conscious
Company in Boulder, Colorado
USA

www.InDevR.com