



VaxArray Influenza Seasonal Hemagglutinin Label Selection Guide

For detailed product information and full assay procedure, please refer to the most recent revision of the VaxArray® Influenza Seasonal Hemagglutinin Operation Manual.

SELECTION GUIDE

Multiple detection labels are available for use with the VaxArray Influenza Seasonal Hemagglutinin Potency Assay (VXI-7100). This guide details the label recommendations for the specific target antigen(s) to be evaluated. For target antigens for which multiple labels are recommended or for unique target antigens other than those listed below, please contact your InDevR Sales Representative to discuss the optimal choice for your testing needs.

Sample Type	VaxArray Influenza Labels								
	A&B pAb Label VXI-7601	A H1-02 mAb Label VXI-7602	A H1-09 mAb Label VXI-7609	A HA-04 mAb Label VXI-7604	A H3-03 mAb Label VXI-7603	A H3-08 mAb Label VXI-7608	A H3-10 mAb Label VXI-7610	B Yam-07 mAb Label VXI-7607	B Vic-06 mAb Label VXI-7606
HA from split A/H1 virus	✓			✓					
HA from split A/H3 virus	✓					✓	✓		
HA from split B/Y-like virus	✓							✓	
HA from split B/V-like virus	✓								✓
HA from A/H1 VLP	✓		✓	✓					
HA from A/H3 VLP	✓					✓	✓		
HA from B/Y-like VLP	✓							✓	
HA from B/V-like VLP	✓								✓
Recombinant A/H1 HA		✓	✓	✓					
Recombinant A/H3 HA					✓		✓		
Recombinant B/Y-like HA	✓							✓	
Recombinant B/V-like HA	✓								✓

Note: Monoclonal labels (those labeled with ‘mAb’ in the table above) are recommended for the most sensitive quantification of protein stability.

For technical support, please contact us at techsupport@indevr.com or 1-303-402-9100.

For sales or quote requests, please contact us at sales@indevr.com or 1-303-402-9100.