



VaxArray Influenza Pandemic Hemagglutinin Label Selection Guide

For detailed product information and full assay procedure, please refer to the most recent revision of the VaxArray® Influenza Pandemic Hemagglutinin Potency Assay Operation Manual.

SELECTION GUIDE

Multiple detection labels are available for use with the VaxArray Influenza Pandemic Hemagglutinin Potency Assay Kit (VXI-7200). This guide details the label recommendations for the specific target antigen(s) to be evaluated. For target antigens for which multiple labels are recommended or for unique target antigens other than those listed below, please contact your InDevR Sales Representative to discuss the optimal choice for your testing needs.

Sample Type	VaxArray Influenza Labels		
	A&B	A HA-04	A H7-05
	pAb Label VXI-7601	mAb Label VXI-7604	mAb Label VXI-7605
HA from split A/H5 virus	✓		
HA from split A/H7 virus	✓		
HA from split A/H9 virus	✓		
HA from A/H5 VLP	✓	✓	
HA from A/H7 VLP	✓		✓
HA from A/H9 VLP	✓	✓	
Recombinant A/H5 HA		✓	
Recombinant A/H7 HA			✓
Recombinant A/H9 HA		✓	

For technical support, please contact us at techsupport@indevr.com or 1-303-402-9100.

For sales or quote requests, please contact us at sales@indevr.com or 1-303-402-9100.