



VaxArray Influenza Seasonal Neuraminidase Label Selection Guide

For detailed product information and full assay procedure, please refer to the most recent revision of the VaxArray® Influenza Potency Assay Operation Manual.

SELECTION GUIDE

Multiple detection labels are available for use with the VaxArray Influenza Seasonal Neuraminidase (NA) Potency Assay. This guide details the label recommendations for the specific target antigen(s) to be evaluated. For target antigens for which multiple labels are recommended or for unique target antigens other than those listed below, please contact your InDevR Sales Representative to discuss the optimal choice for your testing needs.

Sample Type	VaxArray Influenza Labels		
	NA A&B pAb Label VXI-7616	A N1-15 pAb Label VXI-7615	A N2-12 mAb Label VXI-7612
	NA from split A/N1 virus	✓	✓
NA from split A/N2 virus	✓		✓
NA from split B/Y-like virus	✓		
NA from split B/V-like virus	✓		
Recombinant A/N1 NA		✓	
Recombinant A/N2 NA			✓

For technical support, please contact us at techsupport@indevr.com or 1-303-402-9100.

For sales or quote requests, please contact us at sales@indevr.com or 1-303-402-9100.