



### VaxArray Influenza Pandemic Hemagglutinin Label Selection Guide

For detailed product information and full assay procedure, please refer to the most recent revision of the VaxArray® Influenza Potency Assay Operation Manual.

#### SELECTION GUIDE

Multiple detection labels are available for use with the VaxArray Influenza Pandemic Hemagglutinin Potency Assay. This guide details the label recommendations for the specific target antigen(s) to be evaluated. For target antigens for which multiple labels are recommended or for unique target antigens other than those listed below, please contact your InDevR Sales Representative to discuss the optimal choice for your testing needs.

Sample Type	VaxArray Influenza Labels		
	A&B pAb Label VXI-7601	A HA-04 mAb Label VXI-7604	A H7-05 mAb Label VXI-7605
HA from split A/H5 virus	✓		
HA from split A/H7 virus	✓		
HA from split A/H9 virus	✓		
HA from A/H5 VLP	✓	✓	
HA from A/H7 VLP	✓		✓
HA from A/H9 VLP	✓	✓	
Recombinant A/H5 HA		✓	
Recombinant A/H7 HA			✓
Recombinant A/H9 HA		✓	

For technical support, please contact us at [techsupport@indevr.com](mailto:techsupport@indevr.com) or 1-303-402-9100.

For sales or quote requests, please contact us at [sales@indevr.com](mailto:sales@indevr.com) or 1-303-402-9100.